



## Aerotest Simultan Alpha

**For monitoring breathing air in the low pressure range.**

The quality of the breathing air can be tested in accordance to the valid purification standards by means of quantitative measurements of the contaminations. Different regulations can be fulfilled (e.g. DIN EN 12021).

### Quick-check

the measurement time in 5 minutes

### with high accuracy

The measuring device, barbed nipple, can be connected to the low pressure compressed air supply systems.

All components of the Aerotest Simultan Alpha are arranged in a carrying case and ready for use.

Dräger Safety AG & Co. KGaA  
BreathingGas Systems  
Revalstr. 1  
D-23560 Lübeck  
Phone +49 451 882-46 73  
Fax +49 451 882-38 67  
Internet: <http://www.draeger.com>

## Technical Data

Carrying case	Lenght 300 mm; Height 360 mm; Depth 80 mm
Weight	approx 2 kg
Supply pressure	max. 12 bar
Connection	barbed nipple
Flow	0.2 L/min and 4.0 L/min

## Part Number

Aerotest Simultan Alpha, complete	65 27 150
<ul style="list-style-type: none"> <li>- one reducer with hand operation</li> <li>- one measuring device for 4 tubes</li> <li>- one timer</li> <li>- one break-off facility for tubes</li> <li>- one instruction for use manual</li> <li>- one box of 10 Dräger tubes<sup>®</sup>, carbon 100/a-P</li> <li>- one box of 10 Dräger tubes<sup>®</sup>, carbon monoxide 5/a-P</li> <li>- one box of 10 Dräger tubes<sup>®</sup>, oil 10/a-P</li> <li>- one box of 10 Dräger tubes<sup>®</sup>, water vapour 20/a-P</li> <li>- one transport case</li> </ul>	
Aerotest Simultan Alpha - without Dräger tubes <sup>®</sup> ,	65 27 149
<ul style="list-style-type: none"> <li>- one reducer with hand operation</li> <li>- one measuring device for 4 tubes</li> <li>- one timer</li> <li>- one break-off facility for tubes</li> <li>- one instruction for use manual</li> <li>- one transport case</li> </ul>	
Box of 10 Dräger tubes <sup>®</sup> , carbon 100/a-P	67 28 521
Box of 10 Dräger tubes <sup>®</sup> , carbon monoxide 5/a-P	67 28 511
Box of 10 Dräger tubes <sup>®</sup> , oil 10/a-P	67 28 371
Box of 10 Dräger tubes <sup>®</sup> , water vapour 20/a-P	81 03 061